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DESCRIPTION:

NOVATHANE orthoses are made from elastic material with a zipper and zip-puller for easy fit. The materials carry the Standard 100 by Oeko-Tex® label, an independent harmful substances testing and certification system for textile products. NOVATHANE orthoses are single-patient devices. They come in standard sizes and three different positions. The pouring of a semi-rigid polyurethane foam into the orthosis by the health professional permits to fit the patient morphology and to fix the most suitable version.

INDICATIONS:

	NOVARHIZ	Wrist and thumb orthosis Stabilization of the wrist and thumb column.	Rheumatologic and/or neurological disorders as: Rhizarthrosis / Early osteoarthritis / Rheumatoid arthritis; Carpal tunnel syndrome / Wrist tendinitis and tenosynovitis; De Quervain's tenosynovitis; Hypermobile Ehlers Danlos syndrome / Hypermobility spectrum disorders.
	NOVAWRIST	Wrist orthosis with thumb opening Stabilization of the wrist	Rheumatologic and/or neurological disorders as: Early osteoarthritis / Rheumatoid arthritis; Carpal tunnel syndrome / Wrist tendinitis and tenosynovitis; Hypermobile Ehlers Danlos syndrome / Hypermobility spectrum disorders.
	NOVAHAND	Wrist and hand orthosis Hand at rest	Afteintes rhumatologiques et/ou neurologiques suivantes : Flexion contracture / Retraction of the flexor muscles / Hemiplegia; Early osteoarthritis / Rheumatoid arthritis; Carpal tunnel syndrome / Wrist and/or thumb tendinitis and tenosynovitis; Hypermobile Ehlers Danlos syndrome / Hypermobility spectrum disorders.

CLINICAL BENEFIT:

The expected benefit is pain and stiffness reduction, and making certain gestures easier.

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INSTRUCTIONS:

Mixing the polyurethane resin (NOVATHANE A) with the hardener (NOVATHANE B) results in the formation of a microcellular foam designed for manual casting. Both components are packed in two compartments of an airtight bag separated by a plastic ring. When the ring is removed the two components mix together. This liquid mixture is then poured into the opening of the fabric and is evenly spread in the orthosis with a roller and/or by hand (see picture). The orthosis is finally fitted on the patient. The patient maintains the indicated position. The resin becomes semi-rigid in a few minutes.

To illustrate these instructions, a tutorial is available at: <https://www.novatex-medical.com/ortheses-dimmobilisation/>
Regulation (EU) 2020/1149 requires that any worker handling diisocyanates must be trained. The resin containing diisocyanates, it is necessary that the user takes the e-learning course at the following link: <https://www.safeusediisocyanates.eu/>

PRECAUTIONS FOR USE:

- Do not use if the orthosis shows signs of wear;
- Do not modify the orthosis;
- Avoid placing the orthosis directly on a heat source or exposing it to direct sunlight;
- Comply with the medical prescription in terms of treatment duration and indication;
- When applying the orthosis, the skin should be dry, clean and free of cosmetics;
- Use only Novatex Medical resins (NOVATHANE A + NOVATHANE B).

CONTRAINDICATIONS:

- Do not use if the patient is allergic to any of the ingredients.
- Do not apply the orthosis on an open wound or on damaged skin.

UNDESIRABLE EFFECTS:

Skin reactions (redness, itching, skin ulceration, edema, etc.) may occur when using the orthosis. It is recommended to consult your physician in these specific cases. Consult your physician without delay if symptoms worsen or in case of abnormal changes. Any serious incident related to this orthosis must be notified to the manufacturer and the competent authority of the country in which the patient/user is resident.

WARNINGS:

The user must respect the application times of the resin. Do not inhale the resin and avoid contact with the skin and eyes by wearing safety glasses and gloves during use. In case of contact with the skin, rinse with soap and water. In case of contact with the eyes, rinse with water for several minutes. Do not apply if the resin has pierced the fabric.
Do not use the resin if the pouch is damaged. Return any damaged devices to NOVATEX MEDICAL®.
The product must be fitted on the patient by a trained health professional (pharmacist, orthopedist, orthotist, orthopaedic technician, etc.).
Driving vehicles using these orthoses is not recommended.
Do not use the product near an open flame or on a hot surface. This can cause serious health damages.

STORAGE:

Recommended storage conditions for the packed polyurethane system: temperature between 15°C and 30°C. Do not use the product after the expiry date.

CARE:

It is recommended to wash the orthosis by hand with warm water without twisting. The use of softener, stain remover or any chlorinated product is prohibited. Dry in air away from any heat source. Do not tumble dry. Do not iron.

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