DESCRIPTION :

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DESCRIPTION :

NOVASPINE immobilizing braces are made from stretch fabric and feature adjustable front fastening for easy donning. The hook and loop patches to be positioned on the waist belt allow for precise adjustment of the brace. The fabric strips sewn onto the upper part of the brace provide scapular anchorage during moulding, making it easier for the healthcare professional to fit the spinal brace. The materials used are certified according to the Oeko-tex® Standard 100, an independent health label that tests textiles for harmful substances. NOVASPINE braces are single-patient devices. They come in standard sizes and can be made in three lengths corresponding to the apices of the spinous processes, starting from T12, T9 or T6. The health professional casts a rigid or semi-rigid polyurethane foam into the orthosis, which is then moulded to the patient's anatomy in the most appropriate posture.

INDICATIONS :

The NOVASPINE spinal brace is indicated for the post-operative treatment of lumbar disc herniation and for support in patients with lumbar disc disease.

CLINICAL BENEFIT :

The expected clinical benefit is a reduction in pain through limitation of anteroposterior and lateral flexion and as a result of the compressive effect of the brace. Used postoperatively, the motion restriction can also promote wound healing.

COMPOSITION :

Main fabric: 86% Polyamide / 14% Elastane (blue) or 80% Polyamide / 20% Elastane (black); Waist belt: 42% Polyamide / 40% Polyester / 14% Elastane / 4% Polyurethane ; Core: Semi-rigid or rigid polyurethane foam

CONTRAINDICATIONS :

- Do not place the brace in direct contact with an open wound or broken skin (place a suitable sterile interface such as a dressing, gauze or other between the brace and the wound).
- Do not use the brace in case of known allergy to any of the components.
- Do not use on patients for whom compression of the abdomen is contraindicated, in particular pregnant women and patients with circulatory, pulmonary, cardiovascular or skeletal problems.
- Seek medical advice before wearing the brace if you have pelvic floor weakness.

PRECAUTIONS FOR USE :

- Check the brace for signs of wear or deterioration of the components or of the stitching before each use. Do not use if signs of wear are found.
- Do not leave the brace in direct sunlight or exposed to any other source of heat.
- Secure the hook-and-loop fasteners when removing the brace.
- Wear a undershirt under the brace, preferably a seamless one. The brace should be worn between the underwear and the outerwear.
- Follow the physician's instructions regarding the duration of treatment and indication.
- Do not (re)use on other patients.
- Fitting must be carried out by an orthopaedist/orthotist, pharmacist or orthopaedic technician.
- Do not make any modifications to the brace.
- Do not allow the brace to come into contact with cosmetics or harsh products.
- To guarantee traceability, do not remove the label from the brace.
- Do not use the product near an open flame or a hot surface. As the brace is made of synthetic materials, it could be damaged and cause burns or even breathing difficulties if the fumes are inhaled.

UNDESIRABLE EFFECTS :

All braces can exert excessive local pressure if worn too tight. It is therefore essential that such devices are properly fitted and worn. Seek immediate medical advice if symptoms worsen or if any other anomalies occur. Skin reactions (redness, itching, sores, oedema, etc.) may occur during use of the brace. Seek immediate medical advice if this happens. Skin reactions may be due to erosion of the scar (which will remain relatively sensitive for a period of approximately 3 months), skin allergies and intolerances (to a component of the brace or unsuitable detergent), poor fitting of the brace or compression (causing oedema, friction, etc.), or maceration associated with the use of creams and lotions. Any serious incident occurring in connection with the brace must be reported to the manufacturer and to the competent authority in the country where the patient/user is based.

FITTING :

The first time, the brace must be fitted by a trained healthcare professional. This is also an opportunity to train the patient for subsequent use.

CARE AND MAINTENANCE :

To clean the brace, it is recommended to wash it by hand in lukewarm water without twisting and after joining the hook-and-loop fasteners to avoid damaging the fabric. Do not use fabric softener, stain remover or chlorinated products. Allow the brace to air-dry naturally, away from any source of heat. Do not tumble dry or iron. To maintain the grip of the hook-and-loop fasteners, remove any loose threads and lint regularly.

SECTION INTENDED EXCLUSIVELY FOR HEALTHCARE PROFESSIONALS :

FITTING PROCEDURE :

Place the brace on the patient to draw the line along which the back of the brace will be cut as prescribed. Remove the brace from the patient and mix the resin.

Combination of the polyurethane resin (NOVATHANE A or AR) with the hardener (NOVATHANE B or BR) creates a micro-cellular foam designed for the manufacture of hand-cast parts. The two components are packaged in a hermetically sealed pouch and separated by a plastic ring. The two are mixed together once the ring is removed. The liquid mixture is then poured into the opening in the fabric and evenly distributed throughout the brace using the applicator roller and/or by hand (see diagram). Lastly, the brace is placed on the patient. The patient must hold the specified position. After a few minutes, the resin becomes semi-rigid or rigid. The healthcare professional can then trim the back of the brace along the previously drawn line.

A tutorial illustrating these instructions is available at: <https://www.novatex-medical.com/corsets-dimmobilisation/>. According to EU Regulation 2020/1149, appropriate training is a basic necessity for every worker handling diisocyanates. As the resin contains diisocyanates, users must take the e-learning course available at the following link: <https://www.safeusediisocyanates.eu/>

WARNINGS :

- Use only NOVATEX MEDICAL resins (NOVATHANE A + NOVATHANE B or NOVATHANE AR + NOVATHANE BR).
- Comply with the recommended application times for thermosetting resins.
- Do not inhale the resin fumes and avoid all contact with skin and eyes. Users should therefore wear protective goggles and gloves when using the product. In case of contact with the skin, wash with soap and water; in case of contact with eyes, rinse with water for several minutes.
- Do not place the brace on the patient if the resin has flowed through the fabric.
- Do not use the device if it is damaged. Any damaged device should be returned to NOVATEX MEDICAL.
- The product must be fitted by a healthcare professional (pharmacist, orthopaedist/orthotist, orthopaedic technician, etc.).
- During fitting, the skin must be dry, clean and free of any cosmetics.

STORAGE :

Recommended storage conditions for the packaged polyurethane system: temperature between 15°C and 30°C. Do not use the product after the expiry date.

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