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## INSTRUCTIONS FOR USE PRESSURE GARMENTS FOR THE TREATMENT OF SEVERE BURN VICTIMS

### **DESCRIPTION:**

NOVATEX MEDICAL® pressure garments are medical devices custom-made\* in accordance to your individual conditions and needs. They are for use with a single patient in compliance with a medical prescription. The materials carry the Standard 100 by Oeko-Tex® label, an independent harmful substance testing and certification system for textile products.

### **INDICATIONS:**

NOVATEX MEDICAL® pressure garments are indicated in the treatment of extensive, deep second- and third-degree burns when the total surface area exceeds 10% of the body surface, or, when the site of the lesion is on the hands, face, neck or flexion creases.

### **CLINICAL BENEFIT:**

The expected benefits are the improvement of functional signs through pain relief and remobilization, as well as the prevention of scar sequelae through the reduction of scar thickness.

### **COMPOSITION:**

Elastane / Polyamide

### **CONTRAINDICATIONS:**

- Do not apply pressure garments directly on an open wound or damaged skin (unhealed skin);
- Do not apply pressure garments in the case of a known allergy to one of their components.

### **PRECAUTIONS FOR USE:**

- Do not use pressure garments if they show traces of wear and tear;
- The pressure garment's seams must be placed on the outside when putting it on, and the attachment systems must not come into contact with the skin;
- Do not perform any modifications on pressure garments;
- Do not expose yourself to the sun wearing pressure garments;
- Avoid placing the orthosis directly on a heat source or exposing it to direct sunlight;
- Comply with the medical prescription in terms of the length of treatment, the body part concerned, and the indication;
- Do not (re)use on other patients;
- Do not put the orthosis in contact with cosmetics or aggressive products;
- Inform your orthopaedist/orthotist if you are pregnant so they can adapt the pressure around your abdomen;
- If your pressure garment includes a hook and loop fastener, take particular care when putting it on to avoid tearing;
- Do not remove the label on pressure garment for traceability purposes;
- Do not use the product near an open flame or on a hot surface. Since the orthosis consists of synthetic fibres, it can degrade and cause burns.

### **ADVERSE EFFECTS:**

You may feel discomfort when putting the pressure garment on, this is due to the pressure from the device. Consult your doctor without delay if symptoms worsen or in case of unusual issues. Skin reactions (redness, itching, skin ulceration, oedema, etc.) may appear when using a pressure garment. It is recommended to consult your physician in these precise cases. They may be due to scar erosion (the scar remains more fragile for about 3 months), skin allergies and intolerances (to a component or to inappropriate washing products), poor patient matching of garments or pressure (oedema, friction, etc.), a maceration potentially related to cream use or silicone face if applicable. Any serious incident related to this orthosis must be notified to the manufacturer and the competent authority of the country in which the patient/user is resident.

### **USE:**

Before putting it on, the garment can be stretched a bit to make it easier to use. Make sure your skin is clean and dry, with no cosmetics. To avoid damage to the garments, it is advisable to remove jewelry when putting it on. Upon receipt of any new set, the garments should be put on the patient by a healthcare professional (pharmacist, orthopedist, orthotist, orthopedic technician, etc.) for a correct fit to the patient. It is also an opportunity to train the patient for subsequent implementations.

### **CARE:**

For your garment's upkeep, it is recommended to wash it by hand without twisting or in a machine at 30°C. Never use fabric softener, stain remover or any chlorinated products. It should be dried naturally, lying flat, away from sources of heat. Clothes driers and irons should never be used. When washing, fasten the hook and loop systems to avoid damaging the fabric. In order to keep the fasten characteristics of the hook-and-loop devices, maintain them free of threads and fluffs.

### **WARRANTY:**

The warranty (products and labour) covering production, finishings and quality is for a period of six months from the date of delivery. This warranty is valid only in the case of normal use, without damage resulting from external causes and in compliance with the conditions of use and maintenance.

*\* The regulation 2017/745 states that the manufacturer does not affix the CE marking to medical devices if they are custom-made, although they must meet the relevant requirements (article 21).*

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